

## Double mixing Stopped-flow using interrupted flow mode

### I – Introduction

SFM-3000 and SFM-4000 are the ideal systems for double mixing experiments taking full advantage of the precision of independent stepping motors. User can design easily a series of shots at different ageing times so reaction intermediates not visible with single mixing methods become observable.




Fig. 1 : Double mixing set-up.

In a double mixing experiment, two reactants are mixed in mixer 1 and are allowed to age during an user-defined ageing time. After that content of the delay line is mixed with a third reactant so this second reaction can be observed using optical method. Depending on ageing time the mixing is done in continuous flow mode or in interrupted flow mode to cover ageing time range from 2 ms to several minutes.

### II – Experiment

To illustrate double mixing the standard DCIP/ascorbic acid reaction can be used. In mixer 1 dichloro-indophenol (DCIP) is mixed with a low concentration of ascorbic acid (0.5 mM), the mixing is allowed to age for a define time in a 120  $\mu$ l delay line before it is mixed with a higher concentration of ascorbic acid (5 mM). This second reaction is followed in absorbance mode using a TC-100/10 cuvette (1 cm light path).




Fig. 2 : Transients obtained varying ageing time from 55ms to 600 ms.

Total flow rate is kept constant for all ageing times. Transients are recorded with ageing times from 55 ms to 600 ms as shown in figure 2. Each transient is fitted with a single exponential so its amplitude is measured and user can check observed rate constant is constant. As dead time is constant plotting the amplitude of transient versus time will show the first reaction and one can define the rate constant of first reaction

