

Impedance IIIa

Study of the insertion reaction by impedance

Application to the characterization of batteries

Understand to which mechanisms correspond the impedance graphs obtained on batteries.

- 1. Introduction : the Li battery**
- 2. The insertion reaction with restricted diffusion**
 1. Mechanism
 2. Expression of the Faradaic impedance
 3. Equivalent circuit
 4. Experimental data
 5. Other mechanisms
- 3. The insertion reaction with semi-infinite diffusion**
- 4. The insertion reaction with bounded diffusion**

1. Introduction : the Li battery


2. The insertion reaction with restricted diffusion

1. Mechanism
2. Expression of the Faradaic impedance
3. Equivalent circuit
4. Experimental data
5. Other mechanisms

3. The insertion reaction with semi-infinite diffusion


4. The insertion reaction with bounded diffusion

1. Initial state (charged) : $E_{\text{cell}} = E_{\text{OC}}$




E_{cell} = Cell voltage = potential difference between the positive electrode and the negative electrode.

2. Discharge: $E_{\text{cell}} < E_{\text{OC}}$ (charged) (Spontaneous reactions)




3. Charge : $E_{\text{cell}} > E_{\text{OC}}$ (discharged) (Forced reactions)




1. Introduction : the Li battery
2. **The insertion reaction with restricted diffusion**
 1. Mechanism
 2. Expression of the Faradaic impedance
 3. Equivalent circuit
 4. Experimental data
 5. Other mechanisms
3. The insertion reaction with semi-infinite diffusion
4. The insertion reaction with bounded diffusion

Let us consider one electrode of the battery.



No inserted species can flow in the substrate.

At the host/substrate interface, $J_{\langle \text{M} \rangle}(L, t) = 0$



Insertion reaction rate $v_i(t)$

$$v_i(t) = K_i(t) M^+(0, t) \langle \rangle (0, t)$$

with $K_i(t) = k_i \exp[-\alpha_i f E(t)]$

Desinsertion reaction rate $v_d(t)$

$$v_d(t) = K_d(t) \langle M \rangle (0, t)$$

with $K_d(t) = k_d \exp[\alpha_d f E(t)]$

$$f = F/(RT)$$

$$\text{Symmetry factors : } \alpha_i + \alpha_d = 1$$


Global reaction rate $v(t)$

$$v(t) = v_i(t) - v_d(t)$$

$$v(t) = K_i(t) M^+(0, t) \langle \rangle (0, t) - K_d(t) \langle M \rangle (0, t)$$

$$i = i_f(t) + i_c(t) = -Fv(t) + C_{dc} dE/dt$$

Non-steady state regime




$$\frac{\partial \langle M \rangle(x,t)}{\partial t} = D_{\langle M \rangle} \frac{\partial^2 \langle M \rangle(x,t)}{\partial x^2}$$

$$J_{\langle M \rangle}(0,t) = v(t) = -i_f(t)/F$$

$$J_{\langle M \rangle}(L,t) = 0$$

Steady-state regime



$$D_{\langle M \rangle} d^2 \langle M \rangle(x) / dx^2 = 0 \Rightarrow \langle M \rangle(x) = \langle M \rangle$$

In the steady-state regime, the applied potential E is an equilibrium potential $i_f = -Fv = 0$

1. Introduction : the Li battery
2. **The insertion reaction with restricted diffusion**
 1. Mechanism
 2. **Expression of the Faradaic impedance**
 3. Equivalent circuit
 4. Experimental data
 5. Other mechanisms
3. The insertion reaction with semi-infinite diffusion
4. The insertion reaction with bounded diffusion

The Faradaic current is expressed as : $i_f(t) = -Fv(t)$

The Taylor development of the Faradaic current followed by its Laplace transformation leads to an expression of the Faradaic impedance :

$$Z_f(p) = \frac{\Delta E(p)}{\Delta i_f(p)} = R_{ct} + Z_{<M>(p)}$$

R_{ct} is the charge transfer resistance = $\lim_{\omega \rightarrow 0} Z_f(\omega)$

$Z_{<M>}$ is the impedance related to the concentration of the inserted species <M>.

$Z_{<M>}(p)$ can be written :

$$Z_{<M>}(p) = R_{<M>} \frac{\coth \sqrt{\tau_{d<M>} p}}{\sqrt{\tau_{d<M>} p}}$$

$$Z_f(p) = R_{ct} + R_{<M>} \frac{\coth \sqrt{\tau_{d<M>} p}}{\sqrt{\tau_{d<M>} p}}$$

This impedance is equivalent to an electrical circuit.
The impedance of such a circuit can be displayed in a Nyquist diagram.

In this case, $p = j\omega = j2\pi f$ with f the frequency (Hz).

$$\tau_d = L^2/D$$

1. Introduction : the Li battery

2. The insertion reaction with restricted diffusion

1. Mechanism
2. Expression of the Faradaic impedance
3. Equivalent circuit
4. Experimental data
5. Other mechanisms


3. The insertion reaction with semi-infinite diffusion


4. The insertion reaction with bounded diffusion



Faradaic impedance


Electrode impedance





Total measured impedance

Simulated impedance graph using ZSim



Inductance component L (due to cell connections, wires...)

1. Introduction : the Li battery

2. The insertion reaction with restricted diffusion

1. Mechanism
2. Expression of the Faradaic impedance
3. Equivalent circuit
4. Experimental data
5. Other mechanisms

3. The insertion reaction with semi-infinite diffusion

4. The insertion reaction with bounded diffusion

H₂ insertion




Fig. 3 : Experimental impedance diagram for a Pd electrode 100 μ m thick ; potential 115 mV (RHE) ; solution 1 M H₂SO₄ ; temperature 20 °C.

From J.S. Chen, J. -P. Diard, R. Durand, C. Montella, Electrochemistry and Materials Science of Hydrogen Absorption and Adsorption, Electrochemistry Society Meeting, B.E. Conway, G. Jerkiewicz (ed.), San Francisco, (1994) p. 207

H₂ insertion

Obtained on a Pd electrode with $P_{H_2} = 12$ mbar








Fig. 6. Electrochemical impedance diagrams measured on Pd-H at 298 K, $H/Pd = 0.01$ and $E = +60$ mV NHE. (○) Experimental harmonic; (+) experimental non-harmonic with (a) $\Delta t = 50 \mu s$; (b) $\Delta t = 25 \mu s$; (c) $\Delta t = 2.5 \mu s$; (—) model from Eq. (5).

From P. Millet, R. Ngameni, *Electrochim. Acta* 56 (2011) p. 7907

LiFePO₄ battery : the use of ZFit



LiFePO₄ battery : the use of ZFit



Which EC should be chosen ?

The insertion reaction with restricted diffusion

1. Introduction : the Li battery

2. The insertion reaction with restricted diffusion

1. Mechanism
2. Expression of the Faradaic impedance
3. Equivalent circuit
4. Experimental data
5. Other mechanisms

3. The insertion reaction with semi-infinite diffusion

4. The insertion reaction with bounded diffusion

5. The Solid Electrolyte Interphase


Other more complicated insertion reactions are possible

1. Insertion with preliminary electro sorption



Faradaic impedance $Z_f(\omega) = R_{ct} + Z_\theta(\omega)$.

$Z_\theta(\omega)$ is the impedance related to the electroadsorbed species.




Electrode Impedance


2. Insertion with preliminary adsorption



$$\text{Faradaic impedance} = Z_f(\omega) = R_{ct} + Z_{ad}(\omega) + Z_{<M>}(\omega)$$




Simulated Nyquist Plot



Experimental Data


From M.D. Levi et al, J. Electroanal. Chem, 569, 2004



(b)

Mg-ion insertion into the Chevrel electrode ($M_xMo_6S_8$, $0 < x < 2$ for Mg and $0 < x < 4$ for Li)

1. Introduction : the Li battery
2. The insertion reaction with restricted diffusion
 1. Mechanism
 2. Expression of the Faradaic impedance
 3. Equivalent circuit
 4. Experimental data
 5. Other mechanisms
3. The insertion reaction with semi-infinite diffusion
4. The insertion reaction with bounded diffusion



The case for which the length of the host material is very large is a limit case of the direct insertion with restricted diffusion where either L is very large or $D_{\langle M \rangle}$ is very small.


The impedance component of a diffusion in a semi-infinite media is called a Warburg component. It is symbolized by the letter W.

$$Z_w(\omega) = \sigma(1-j)/(j\omega)^{1/2}$$


$$Z_f(\omega) = R_{ct} + Z_w(j\omega)$$

Electrode impedance

Such a circuit is called a **Randles** circuit




The vertical component on the Nyquist diagram is shifted to the much lower Frequencies, the time constant τ_{d1} ($= L^2/D_{<M>}$) is much larger.



$R_{d1} = 0,2418 \Omega$ instead of $0,02418 \Omega$, $\tau_{d1} = 100,17$ s instead of $67,17$ s

Simulated impedance graph using ZSim



Randles circuit + electrolyte resistance


1. Introduction : the Li battery
2. The insertion reaction with restricted diffusion
 1. Mechanism
 2. Expression of the Faradaic impedance
 3. Equivalent circuit
 4. Experimental data
 5. Other mechanisms
3. The insertion reaction with semi-infinite diffusion
4. The insertion reaction with bounded diffusion

The impedance component of a bounded diffusion is a W_δ component.

$$Z_{w\delta}(j\omega) = R_d \operatorname{th}(\tau_d j\omega)^{1/2} / (\tau_d j\omega)^{1/2}$$


Electrode impedance

$$Z_f(j\omega) = R_{ct} + Z_{w\delta}(j\omega)$$



Faradaic impedance

Simulated impedance graph using ZSim



- The reaction that takes place in a battery undergoing a charge or discharge is an insertion.
- The insertion involves a diffusion of the inserted species in three different conditions : restricted, semi-infinite, bounded.
- The insertion mechanism can be direct or involve a preliminary electrosorption or adsorption.
- For all these conditions, we now know :
 - The expression of the Faradaic impedance
 - The equivalent circuit
 - The Nyquist diagrams of the impedance

- Having this knowledge, we now know how to interpret impedance data obtained on batteries.
- The next tutorial Impedance IIIb will show what useful characteristics of the batteries can be obtained from the impedance data.

Bio-Logic Website : www.bio-logic.info

<http://www.bio-logic.info/potentiostat/notesifil.html>

Faradaic Impedances

Mathematica Player files :

Direct Insertion

<http://www.bio-logic.info/potentiostat/notes/20080131%20-%20Insertion1-ZmmaP.nbp>

Electrosorption + Insertion

<http://www.bio-logic.info/potentiostat/notes/20080307%20-%20IndirectInsertion1-mmaP.nbp>

Adsorption + Insertion

<http://www.bio-logic.info/potentiostat/notes/20080312%20-%20IndirectInsertion2-mmaP.nbp>

Equivalent Circuits

<http://www.bio-logic.info/potentiostat/iecl/20101004%20-%20RandlesCircuit-mmaP.nbp>