

The use of sweep scan in Scanning Kelvin Probe (SKP) experiments

I. Introduction

Scanning Kelvin Probe (SKP) measures the contact potential difference, or Volta potential, between the SKP probe and the sample under study. It is possible to correlate the Volta potential measured with the work function or the corrosion potential of the sample under study.

A limiting factor in SKP, and other scanning probe experiments, can be the length of the experiment. This is particularly the case when large scan areas, of a few centimeters or more, are of interest. In these cases, using step scan, in which the probe stops to measure at each point, can lead to prohibitively long experiments. If the Volta potential is measured continuously, as in sweep scan mode, the SKP experiment duration can be noticeably reduced, making it highly applicable to the measurement of large areas.

This note uses the M470 to perform SKP measurements of a laser marked anodized aluminium sample, provided by Renthal Ltd. Laser marking is an industrially important technique used to inscribe information and decoration on parts. When used on anodized aluminium the laser marker removes the anodized layer exposing the underlying aluminium [1]. The laser marking was done at different intensities and frequencies to form a grid of varying exposures. Measurements were performed in step and sweep scan to demonstrate the applicability of sweep scan to large area SKP experiments.

II. Sweep vs. step scan

Scans for all dc techniques on the M470 can be performed in either step scan or sweep scan mode. In step scan mode the probe measures the first point in each line before moving to the next point, pausing, and performing the next measurement.

This move-pause-measure cycle is repeated at each point until the line is complete, this is illustrated in Fig. 1a. Because the probe pauses at each point multiple samples are measured and averaged. Step scan is used, therefore, for measurements with the lowest noise. In sweep scan, Fig. 1b, the probe does not stop and measures continuously throughout each line. Because the velocity of the probe is accurately known measurements are taken at set time intervals and converted into the position data of interest. As will be shown in this note sweep scan is beneficial when the shortest experiments are needed, for example when the size of the scan area could lead to prohibitively long experiment durations.




Figure 1 : Step scan (a) and sweep scan (b) modes are illustrated.

III. Experimental

Laser marking of the blue anodized aluminium sample was performed at a speed of 500 mm/s. In the x direction the frequency increased from 30 kHz to 100 kHz, while in y the power increased from 5% to 100%. The result is a grid of 15 squares in x, and 17 squares in y, Fig. 2.

The laser marked sample was mounted on the TriCell™. To accommodate the large size of the sample the 32 mm sample holder was removed from the TriCell™ base allowing it to rest directly on the Perspex base. The sample is initially levelled, by placing a spirit level directly on top of it. Once the sample was satisfactorily levelled with the spirit level the Capacitive Height Measurement (CHM) was used to make fine adjustments to the sample level. In the CHM configuration a real time output of the probe to sample distance is displayed. By moving the probe to the four corners of the measurement area it is possible to level the sample by obtaining similar CHM outputs at all locations. A crocodile clip was used to make electrical connection to the sample. A 500 µm tungsten probe (U-SKP370/1) positioned within 100 µm of the sample surface, was used for all measurements.




Figure 2 : The laser marked anodized aluminium sample.

Measurements comparing sweep scan and step scan were performed over a 1 cm x 1 cm area, covering nine different laser marked squares. Both measurements used a 100 µm step size in x and y, and a velocity 500 µm/s.

To illustrate the effect of sweep rate on the measurement quality a second sweep scan experiment was performed over a 1 cm x 1 cm area. Only the sweep velocity was changed to 1000 µm/s.

The signal measured in SKP can be dependent on sample topography. The removal of the anodized layer during laser etching causes the etched areas to have a lower topography than those with the anodizing still intact. To rule out the effect of topography from the etching on the final SKP measurements a 1 cm x 1 cm measurement was performed with height tracking. In this test the sample topography was first measured with CHM, and then input into the SKP measurement to maintain the probe to sample distance. The step size was 100 µm in x and y, and the velocity was 500 µm/s.

Finally, a sweep scan measurement over an area of 7.5 cm x 8.4 cm, to cover the entire laser marked grid, was performed. The step size and velocity were 250 µm, and 500 µm/s respectively.

Final area maps were produced using Gwyddion [2] to rotate them to match the orientation of the laser marked piece.

IV. Results

1. Sweep scan vs step scan

Fig. 3 shows the initial sweep and step scan experiments. These have been rotated 90° counter clockwise from the measurement direction to match the orientation of the laser marked sample. There is little difference between the Volta potentials measured in either mode. In the axis of the sweep scan, however, the boundary between the anodized layer and the laser marked squares are sharper in the step scan than in the sweep scan measurement. This is a result of the data averaging which occurs in the step scan measurement, improving the final signal. The signal improvement, comes at a cost. In this example the 1 cm x 1 cm area scan took 42.5 minutes in sweep scan mode, and 4 hours and 56 minutes in step scan mode, an almost seven-fold increase in experiment duration.




Figure 3 : SKP measurement of the same area of laser marked sample using the (a) sweep scan, and (b) step scan modes.

2. Effect of sweep rate

When using sweep scan it is important to note that while faster sweep rates reduce the experiment duration, they can also exaggerate the effect of noise, and degradation of the measurement signal can occur, as demonstrated in Fig. 4. Fig. 4a is the sweep scan performed at 500 $\mu\text{m/s}$, while Fig. 4b is performed at 1000 $\mu\text{m/s}$. While using the faster sweep velocity has reduced the measurement duration to 22.25 minutes, almost half as long as the sweep scan at 500 $\mu\text{m/s}$, the border between the laser marked square and anodized aluminium is further blurred. Furthermore, more variation in the signal measured within each laser marked square is seen. When using sweep scan, therefore, users must consider how the selected sweep rate will affect the final signal. By using a more conservative sweep rate users can still take advantage of the reduced experiment duration compared to step scan, while achieving sharper, higher quality measurements.




Figure 4 : SKP sweep scan of the same area of a laser marked sample at (a) 500, and (b) 1000 $\mu\text{m/s}$

3. Topography effect

Before measuring the full laser marked grid it was necessary to rule out the possibility that the changes in Volta potential measured are related to changes in topography as a result of the laser marking of each square. To do this the topography of the same 1 cm x 1 cm area was measured using the CHM mode. The resulting CHM map is shown in Fig. 5. This map has been rotated 90° counter clockwise from the measurement direction to match the orientation of the laser marked sample. The CHM map shows the maximum height difference between the anodized and etched regions of approximately 20 μm , with the largest height differences at the squares with the most intense etching. Using this CHM map as a topography input for height tracking, it was then possible to re-measure the sample with SKP, this time with the probe to sample distance maintained throughout. The result of this measurement is compared to the original constant height SKP measurement in Fig. 6. When these two maps are compared no apparent difference in

quality or magnitude of the measurement can be seen, therefore a topography effect can be ruled out.




Figure 5 : CHM topography measurement of 1 cm x 1 cm area of laser marked sample.




Figure 6 : SKP step scan measurement of the same area of a laser marked sample in (a) constant height and (b) height tracking mode.

4. Large area sweep scan

As a demonstration of the applicability of the use of sweep scan the final

measurement was of the entire laser marked grid. Using a 500 μm/s sweep rate, for the highest quality measurement with the shortest experiment duration, the entire 7.5 cm x 8.4 cm grid was measured in just over 15.5 hours. Extrapolating from the original step scan measurement to measure an area of this size, with the same number of points, over two days would be needed, a length which could exclude the use of SKP to measure the sample.

The final sweep scan is shown in Fig. 7. This has been rotated 180° to match the orientation of the sample. The sample has been levelled using CHM as outlined previously, to avoid the effect of sample tilt which can be more prominent over such large areas. In an SKP measurement the measured Volta potential can be directly related to the corrosion potential of a sample [3-5], with lower Volta potentials indicating a lower corrosion potential. Bearing this in mind, therefore, it can be seen that the laser marked regions have a lower corrosion potential than the intact anodized layer, as would be predicted. The corrosion potential also increases as the intensity of the laser marking decreases. Finally, it can be seen that potential changes can be detected in laser marked regions even where a visual change in the anodized layer has not occurred.




Figure 7 : SKP measurement of the full 7.5 cm x 8.4 cm laser marked grid sample.

V. Conclusion

This note demonstrates the applicability of sweep scan mode, particularly for SKP measurements. While measurements with sweep scan can have a lower signal quality than those performed using step scan mode, this can be an acceptable compromise to noticeably reduce the experiment duration. This is particularly the case when large scan areas are of interest.

VI. Acknowledgements

We warmly thank [Renthal Ltd](#) for providing the samples used in this work.

VII. Appendix: Sweep scan settings

Table 1: Signal Conditioning and Backing Potential settings used for 500 $\mu\text{m}/\text{s}$ sweep scan SKP experiment.

Signal Conditioning	
Gain	100
Full Scale Sensitivity	5 mV
Output Time Constant	0.1
Vibration Amplitude	30 μm
Vibration Frequency	80 Hz
Reference Phase	Manually set
Backing Potential	
P-Proportional	0.700
I-Integral	0.005
D-Derivative	0.000
PID Loop Rate (Hz)	25
ADC Sample Rate (Hz)	1000
Samples Averaged per Loop	20
Max PID Step per Loop	0.5
Max Control Output	5.0
Min Control Output	-5.0

When configuring SKP sweep scan experiments the Signal Conditioning and Backing Potential settings need to be updated from those used in step scan mode to ensure data sampling keeps up with the scanning speed. If sampling is too slow

the signal will lag behind the real feature position, and appear unnaturally smooth. If too fast the data can oscillate resulting in excessively noisy data. While the exact settings, particularly the PID values, may vary between laboratories and samples, the Signal Conditioning and Backing Potential settings used for the 500 $\mu\text{m}/\text{s}$ sweep scan are given below as a reference for users.

REFERENCES

1. L. Lazov, H. Deneva, P. Narica, Environment. Technology. Resources. Proceedings of the 10th International Scientific and Practical Conference. Volume I, 2015, 108-115
2. Gwyddion data analysis software <http://gwyddion.net/>
3. M. Stratmann, H. Streckel, Corrosion Science, 1990, 30, 681 – 696
4. M. Stratmann, H. Streckel, Corrosion Science, 1990, 30, 697 – 714
5. M. Stratmann, H. Streckel, K. T. Kim and S. Crockett, Corrosion Science, 1990, 30, 715 – 734

Revised in 06/2021

Scientific articles are regularly added to BioLogic's Learning Center:

<https://www.biologic.net/topics/>

To access all BioLogic's application notes, technical notes and white papers click on Documentation:

<https://www.biologic.net/documents/>